



Monte Carlo for CyberKnife® Radiosurgery with the Incise™ Multileaf Collimator

J R DOOLEY¹, J M NOLL¹, W KILBY¹, W FONG¹, T YEUNG¹, L M GOGGIN¹, D SPELLMAN¹, J S LI², C-M MA², C R MAURER JR¹

¹ – ACCURAY INC., SUNNYVALE, CA, USA

² – FOX-CHASE CANCER CENTER, PHILADELPHIA, PA, USA

DISCLOSURES

- Eight of the authors are employees of Accuray, Inc.
- Portions of this work were completed through a contract between the Fox-Chase Cancer Center and Accuray, Inc.

The views expressed in this presentation are those of the presenters and do not necessarily reflect the views or policies of Accuray Incorporated or its subsidiaries. No official endorsement by Accuray Incorporated or any of its subsidiaries of any vendor, products or services contained in this presentation is intended or should be inferred.

CyberKnife® with InCise™ Multileaf Collimator

- The CyberKnife is an image-guided therapeutic radiation delivery system with a linear accelerator mounted on a robot
- Beams are collimated with fixed cones, the Iris™ variable aperture collimator, or the InCise multileaf collimator (MLC)
- The MLC has a 11.5 x 10.0cm field size with 26 0.385cm leaf pairs
- Full interdigitation and 100% overtravel are provided





Overview

- Monte Carlo dose model
 - Source model and sampling
 - Transport through the MLC
 - Transport through and dose deposition in the patient
- Accuracy testing
 - Planar single beam comparisons with measurement
 - Point dose comparisons with measurements of composed plans
- Implemented dose calculation times


Source model

- A virtual source model with a single source coincident with the Linac target is defined
- Three probability distribution recreate the photon phase space for the linac
 - Source position distribution: a Gaussian with user set full width half-max
 - Fluence distribution:
 - Measure 0, 45, 90, 135 degree open field profiles with MLC removed
 - Planar distribution interpolates profiles in polar coordinates
 - Energy distribution: user selects Geant pre-calculated spectrum based on linac energy rating
- Gaussian FWHM and Energy spectrum settings are optimized iteratively by comparing MC dose calculations with commissioning beam data measurements



Transport through the MLC

- Leaves are flat sided and move on linear trajectories
- To reduce leakage, the collimator bank is 9cm high and is tilted 0.5°
- Leaf tip has a trapezoidal edge
- Particle transmission is checked in 4 planes that define the trapezoid
 - Set transmission probability based on planes intersected by particle trajectory
- Collimator scatter is not modelled




	PLANE				Transmission Probability
	A	B	C	D	
X					0.32
				X	0.32
		X	X		0.01
		X			0.20
			X		0.20
X	X				0.003
X	X	X			0.003
	X	X	X		0.003
		X	X		0.003
					1.0

Transport through the collimator

Beam Hardening

- Energy for photons that intersect the MLC leaves are sampled from a hardened energy spectrum
- The hardened energy spectrum is calculated from the open spectrum, by assuming 3cm of transmission through tungsten



Transport through the patient or phantom


- Particle transport uses pre-simulated electron track data for mono-energetic photons that interact with water
 - Energy bins range from 25 keV to 7.7 MeV
 - 1000 to 10000 photon simulations in each energy bin
 - Track lengths scaled proportionally based on local density
- Variance reduction techniques include:
 - Electron track repeating
 - Forced photon interactions
 - Photon splitting
 - Russian roulette on collimator attenuated and patient scattered photons
 - Energy cutoffs (10 keV for photons, 700 keV for electrons and positrons)
- Ma, et al, “Implementation of Monte Carlo dose calculation for CyberKnife treatment planning”, J. of Physics: Conf. Series 102 (2008)
- Li, et al, “Validation of a Monte Carlo dose calculation tool for radiotherapy”, Phys. Med. Biol. 45: 2969-85 (2000)

Transport through the patient or phantom

Handling electrons and positrons below E_{cut}

- For the original implementation, if the particle energy falls below E_{cut} during the final step of the pre-generated particle track then the remaining energy is deposited in the current voxel
- For the new version, linearly extrapolate the final track step using stopping-power ratio and local density to attenuate the energy
 - Prevents artificial high-dose sparkles in low density

Original



New

Accuracy Tests

Comparisons of calculated and measured dose

- Single beam tests
 - PDD in homogeneous and heterogeneous phantoms
 - Planar analysis of rectangular and irregularly shaped beams in homogeneous and heterogeneous phantoms
- Composed plans in homogeneous and heterogeneous phantoms

Comparison of PDD in water phantoms

PTW60018 diode vs. PTW μ Diamond detector

- Both detectors give good agreement for field sizes greater than 2 x 2 cm
- μ Diamond has better agreement for small field sizes and deeper depths
- Tabulate ratio of measured to calculated dose

PTW60018 Diode


Field Size (cm)	5cm Depth	10cm Depth	20cm Depth	29.5cm Depth
0.76 x 0.77	0.995	0.989	0.969	0.957
1.54 x 1.54	0.999	0.994	0.980	0.973
2.30 x 2.31	0.996	0.998	0.993	0.993

μ Diamond Detector

Field Size (cm)	5cm Depth	10cm Depth	20cm Depth	29.5cm Depth
0.76 x 0.77	1.007	1.005	1.008	0.988
1.54 x 1.54	1.002	1.007	0.999	0.996
2.30 x 2.31	1.008	1.008	1.010	1.004

Dose Accuracy Tests




PDD in a Water / Lung / Water phantom



Dose Accuracy Tests







Planar gamma analysis of film in Water / Lung / Water phantom

- 3 rectangular shapes (0.76x0.77cm, 1.54x1.54cm, and 4.62x4.62cm) and 9 irregular shapes
- Same water / lung / water phantom as the PDD measurement
- Film measurements recorded 9cm (mid-lung), 13.5cm (0.5cm of buildup), and 14.0cm (1.0cm of buildup) deep
- Gamma analysis at 2% dose difference and 0.2cm distance to agreement for all pixels in the plane with dose 50% or more of the maximum dose



Dose Accuracy Tests

Planar gamma analysis of film in Water / Lung / Water phantom



Dose Calculation Tests

Planar gamma analysis (2.0% / 0.2cm) with film measurement

Shape	9 cm	13.5 cm (2.0mm voxel)	13.5 cm (0.5mm voxel)	14.0 cm
0.76 x 0.77 cm	100%	86.7%	97.7%	91.3%
1.54 x 1.54 cm	94.8%	69.9%	95.7%	100.0%
4.62 x 4.62 cm	100%	97.8%		100.0%
HN Ring	98.1%	69.3%	95.7%	100.0%
HN Random	100%	99.3%		99.8%
HN Conf. Avoid	100%	95.8%		100.0%
Lung Ring	100%	71.5%	90.2%	99.1%
Lung Random	100%	53.0%	81.5%	100.0%
Lung Conformal	100%	100.0%		100.0%
Prostate Ring	95.4%	75.2%	94.4%	95.7%
Prostate Random	100%	63.9%	97.0%	98.7%
Prostate Conf. Avoid	100%	75.0%	100.0%	100.0%

Composed Plans

Point dose comparisons of measurement and calculation

- 3 relatively homogeneous phantoms (9 plans)
- 2 anthropomorphic chest heterogeneous phantoms with low density regions (3 plans)


Phantom	Treatment Plan Segment Types	Measure (cGy)	Monte Carlo (cGy)	MC Diff
Head	A 30, Conf	408.3	403.7	-1.1%
Head	B 30, Conf	4085.4	4033.8	-1.3%
Head	C 35, Conf	406.7	402.8	-1.0%
Head	D 114, Many	393.4	389.9	-0.9%
Head	E 88, Many	393.3	393.5	+0.1%
DQAA16	F 13, Conf	414.6	405.7	-2.1%
DQAA16	G 55, Many	400.0	390.4	-2.4%
DQA	H 13, Conf	414.1	408.9	-1.3%
Cheese ++	I 37, Conf	715.1	694.3	-2.9%
Heterogeneous Phantom Measurements				
XLT Lung	J 22, Conf	341.3	341.0	-0.1%
XLT Lung	K 100, Many	338.7	342.4	+1.1%
RPC Lung (*)	L 126, Many	562.3	568	+1.0%
RPC Lung (*)	M 126, Many	550.4	562	+2.1%

(*) TLD measurement, otherwise ion chamber measurement

Comparison with EGSnrc





- Composed plan comparisons performed by Inselspital – University Hospital Bern, using an independently developed MC dose calculation framework (IDC)¹
- 7 lung SBRT cases evaluated
 - Mean dose difference (EGSnrc-Precision)/Precision
 - 2% / 0.1cm gamma comparison for voxels greater than 10.0% of maximum dose

Case	PTV mean dose diff (%)	Lung V20 (%)	3D gamma passing rate (%)
1	1.4	1.5	99.7
2	-0.4	-1.5	97.9
3	0.0	0.0	99.5
4	2.3	0.0	99.1
5	-1.0	0.9	97.6
6	0.4	0.3	99.1
7	0.9	0.0	99.5



1. P.-H. Mackeprang, ... M.K. Fix, P. Manser, (submitted to PMB)

Results



Calculation Time Analysis

- In the Accuray Precision™ Treatment Planning System, Monte Carlo runs in 23 parallel threads on dual 6-core hyper-threaded CPUs
- Calculation times depend on the number of photons sampled
- The number of photons sampled (n_{hist}) is calculated from:
 - The number of voxels
 - The desired calculation uncertainty (u)
 - The mean equivalent square size of the MLC shapes (\overline{EQS})

$$n_{hist} = H_{ref}' \left(\frac{1}{u}\right)^2 \left(\frac{V}{v}\right)^{2/3} \frac{1}{\overline{EQS^2}}$$

- $H_{ref}' = 1.4E6$, $V =$ target volume, $v =$ voxel volume

This algorithm is implemented in Accuray Precision v1.1.1


Sample Calculation Times

- All plans are lung SBRT using actual patient CT and volumes
- All cases calculated at medium resolution (256 x 256 x number of slices) with a requested uncertainty of 2.0%
- Time is from calculation button clicked to display of isodoses and DVH

Case	PTV (cm ³)	Voxel size at medium res (mm)	Number of Beams	Achieved Uncertainty (%)	Calculation Time (s)
1	93	2.2 x 1.0 x 2.2	23	1.9	142
2	14	1.7 x 1.5 x 1.7	24	2.1	61
3	22	2.2 x 1.3 x 2.2	22	1.8	57
4	22	2.3 x 1.5 x 2.3	19	2.2	52
5	30	1.7 x 1.0 x 1.7	21	1.9	101
6	39	2.0 x 2.5 x 2.0	24	1.7	44
7	83	2.0 x 1.3 x 2.0	21	1.9	140

Calculations performed using Precision v1.1.1 on a Dell T7910, 2 x 2.40Gz CPU, 64 GB RAM

Multi-Threading



- Initial version allocated equal numbers of beams to each thread
- Latest version allocates equal number of sampled photons to each thread
- Thread balancing has improved significantly, but further improvement is possible by balancing based on estimated number of transported photons

Calculations performed using Accuray Precision v1.1.0 (left) vs. v1.1.1 (right)

Conclusion

- Most recent version of Accuray Precision treatment planning system includes a Monte Carlo dose calculation algorithm that can calculate dose for MLC collimated beams
- The algorithm computes dose with clinically acceptable agreement with measurement in both homogeneous and heterogeneous media
- The algorithm has a balanced parallel implementation and takes advantage of multiple variance reduction techniques that lead to clinically realizable calculation times